



Single Use Grasping Forceps

PolyGrab_M

Advanced tripod configuration in a convenient single use design





Wide opening maximizes grasping capability in addition to an innovative ringed tip designed to minimize tissue damage.

Product Features

- Optimized for polyp and foreign body retrieval, the pre-sterilized, single use Polygrab has been designed to provide convenience and efficiency.
- The Polygrab's grasping section incorporates OLYMPUS' acclaimed tripod design, with wide opening prongs and ringed distal tips to ensure excellent retrieval capability.
- The ringed distal tips of the Polygrab's tripod prongs are designed to minimize tissue damage.



Specifications

Model	Working length	Opening wigth	Grasper design	Minimum channel diameter
FG-600U	2300 mm	20 mm	Tripod	2.8 mm

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

