


EndoTherapy for small intestinal procedures



Versatile EndoTherapy accessories from OLYMPUS offer exceptional support for small intestinal procedures, advancing enteroscopy into the future

Sampling

Biopsy Forceps

EndoJaw™

Single Use

Employment of an improved sheath section enables smooth insertion into a scope channel. These forceps are especially designed to prevent damage to the interior of a scope channel. Moreover, the double-step blades on the cups ensure sharp cutting quality, and swing jaw mechanism enables easier cup positioning for biopsy sites approached tangentially.



FB-210U	Alligator Jaw without needle
FB-220U	Alligator Jaw with needle
FB-230U	Standard oval without needle
FB-240U	Standard oval with needle
Minimum Channel Size	2.8 mm
Working Length	2300 mm
Quantity	20 pcs / Box

Reusable

A hole of approximately 1 mm diameter at the center of the cups designed to help minimize tissue damage and provide large sample yield.



FB-28Y-1	Standard oval without needle
FB-28Z-1	Standard oval without needle
Minimum Channel Size	2.8 mm
Working Length	FB-28Y-1 : 2700 mm FB-28Z-1 : 3000 mm
Quantity	1 pc / Box

Foreign Body Retrieval

Retrieval Basket

Especially designed for foreign body removal in small intestine, this retrieval basket offers reliable performance in retrieving resected tissue and other foreign objects with 6 wire basket and rotatable function.

Single Use



FG-460YR	6 wire basket
Minimum Channel Size	2.8 mm
Working Length	2700 mm
Basket Design	6 wire basket
Rotatability	Yes
Quantity	1 pc / Box

Polypectomy

Hot Biopsy Forceps

EndoJaw Hot™

Disposable hot biopsy forceps with the slim outer diameter, the proprietary sheath design and improved flexibility ensure smooth insertion into the endoscope channel.

Single Use




FD-210U	Alligator Jaw
FD-230U	Standard oval
Minimum Channel Size	2.8 mm
Working Length	2300 mm
Quantity	5 pcs / Box

Electrosurgical Snares

SnareMaster™

Single Use


Disposable design with integrated handle for ease of use.



SD-221U-25	Crescent type
Minimum Channel Size	2.0 mm
Working Length	2300 mm
Loop Opening Width	25 mm
Quantity	10 pcs / Box

Reusable

Especially designed for small intestinal procedure.



SD-310Z-25	Oval type
Minimum Channel Size	2.8 mm
Working Length	3000 mm
Loop Opening Width	25 mm
Quantity	2 Wires and 2 tube sheaths



With the introduction of an innovative, yet surprisingly simple Single Balloon Enteroscope System, OLYMPUS has literally opened up the small intestine to enteroscopy, helping this promising field take a quantum leap forward. Our sophisticated line of versatile EndoTherapy accessories combine with the advanced capabilities of the Single Balloon Enteroscope to create a synergy that will put you at the cutting edge of enteroscopy.

Hemostasis

Rotatable Clip Fixing Device

QuickClip2™

Especially designed for enteroscopy procedures, the enhanced rotation feature facilitates easy and effective targeting of lesions in the small intestines.

Single Use



HX-201YR-135	
Minimum Channel Size	2.8 mm
Working Length	2700 mm
Clip Jaw Angle	135°
Rotatability	Yes
Clip Opening Width	Standard
Quantity	5 pcs / Box

Bipolar Probes

BiCoag™

The single use BiCOAG hemostasis probe provides reliable hemostasis through bipolar coagulation technology.

Single Use



CD-B620LA	Coaxial
CD-B622LA	Fixed Pin
Minimum Channel Size	2.8 mm
Working Length	3500 mm
Quantity	1 pc / Box

Injection Needle

InjectorForce Max™

Thanks to improved needle extrusion capability, enhanced click response, and excellent liquid feeding, the InjectorForce Max maximizes endoscopic injection performance.

Single Use



NM-400Y-0423	Short bevel
Minimum Channel Size	2.8 mm
Working Length	2700 mm
Needle Length	4 mm
Gauge	23 G
Quantity	5 pcs / Box

Other Accessories

Spray Catheters

Especially designed for spraying dye, these catheters enable circular spraying of the liquid from the distal end.

Reusable



PW-1V-1	Standard type
PW-5V-1	Spray type
Minimum Channel Size	2.8 mm
Working Length	2400 mm
Feature	PW-1V-1 : Inside Chip (Routine washing) PW-5V-1 : Screw tip (Staining)
Quantity	PW-1V-1 : 2 pcs / Box (A), 1 pc / Box (B) PW-5V-1 : 1 pc / Box

Cleaning Brushes

Designed for one-time use, this brush is ideal for cleaning suction channels and openings of compatible OLYMPUS endoscopes.

Single Use



BW-412T	Channel cleaning brush and channel-opening cleaning brush
Working Length	2200 mm
Quantity	50 pcs / Box

Reusable



BW-9Y	Coil sheath channel cleaning brush
MH-507	Channel-port cleaning brush
Working Length	2850 mm (BW-9Y)
Quantity	1 pc / box

Disposable Distal Attachment

When attached to the distal end of a scope, this transparent cap assists in observation of the mucosa and helps achieve a good field of view.

Single Use



D-201-10704	Soft Cap, Straight type
Max. Outer Diameter	11.35 mm
Distal Length from Scope	4 mm
Applicable Scope	SIF-Q180, Q260*
Quantity	10 pcs / box

*The above mentioned endoscopes are not available in some areas

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.



OLYMPUS MEDICAL SYSTEMS CORP.
Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan

For a complete listing of
sales and distribution locations visit:
www.olympus.com